

VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518013



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
5 ANS
5
5 YEARS
GUARANTEE

07/06/2026 17:21:33 – For information purposes only, can be subject to change

VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518013

=====

Description

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in nylon pouch.

Adult (M), cotton, velcro, navy: 518012

Adult (M), cotton, velcro, red: 518013

Adult (M), navy, with straps: 518112

Adult (M), nylon, velcro, grey: 518011

Adult (M), red with straps: 518113

Adult (L), nylon, velcro, grey: 518002

Child/small adult (S), cotton, velcro, motifs : 518022

 07/06/2026 17:21:33 – For information purposes only, can be subject to change

VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518013

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446008888	HS Code	90189010
Warranty (months)	60	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.757
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/06/2026 17:21:33 – For information purposes only, can be subject to change