


# VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE

Ref. 518112



CE  
0459

DISPOSITIF MEDICAL  
CLASSE  
**I**  
m  
MEDICAL DEVICE

GARANTIE  
**5**  
ANS  
5 YEARS  
GUARANTEE



📅 07/06/2026 16:25:45 – For information purposes only, can be subject to change

# VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518112



## Description

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in nylon pouch.

Adult (M), cotton, velcro, navy: 518012

Adult (M), cotton, velcro, red: 518013

Adult (M), navy, with straps: 518112

Adult (M), nylon, velcro, grey: 518011

Adult (M), red with straps: 518113

Adult (L), nylon, velcro, grey: 518002

Child/small adult (S), cotton, velcro, motifs : 518022

# VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518112

## Features

|                                    |               |                           |                               |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                         | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                           | 3700446002077 | HS Code                   | 90189010                      |
| Warranty (months)                  | 60            | Country                   | France                        |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |                               |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc                            |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.813                         |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |                               |
| Minimum temperature                | -20           | Max temperature           | 70                            |
| Minimum humidity                   | 20            | Maximum humidity          | 85                            |
| <b>OTHER FEATURES</b>              |               |                           |                               |
| Packed unit / Quantity             | 1             |                           |                               |

📅 07/06/2026 16:25:45 – For information purposes only, can be subject to change