

Ref. 332039



CE



09/06/2026 00:07:39 – For information purposes only, can be subject to change

Ref. 332039


Description

Cotton cuff with highly resistant rubber bladder.
 Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.
 Washable at 30-40°C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446008789	HS Code	9018901000
Warranty (months)	6	Country	Indonesia

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	15	Packed unit / Height (cm)	2
Packed unit / Length (cm)	13	Packed unit / Weight	0.062

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

 09/06/2026 00:07:39 – For information purposes only, can be subject to change