

**Ref. 332049**



CE



09/06/2026 00:10:17 – For information purposes only, can be subject to change

---

**Ref. 332049**


## Description

Cotton cuff with highly resistant rubber bladder.  
 Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.  
 Washable at 30-40°C.

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446008772	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	2
Packed unit / Length (cm)	13	Packed unit / Weight	0.034
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 09/06/2026 00:10:17 – For information purposes only, can be subject to change