

## NIBP CUFF IMMERSIBLE CLINIC®

Ref. BM103S



 08/06/2026 23:58:09 – For information purposes only, can be subject to change

# NIBP CUFF IMMERSIBLE CLINIC®

Ref. BM103S

The comfortable, waterproof blood pressure cuff!

## Description

The patient-friendly PNI CLINIC® cuff.

Suitable for all types of patients, with 6 sizes available (5 cuffs and one thigh cuff)  
 Made from TPU, a material that feels pleasant against the skin  
 Easy to clean: can be submerged in water

Available in 5 sizes:

- BM101S infant size
- BM102S child size
- BM103S small adult size
- BM104S adult size
- BM105S large adult size

SPENGLER MEDICAL IS COMMITTED TO SUSTAINABILITY!

For CLINIC® cuffs and thigh cuffs, SPENGLER MEDICAL has replaced plastic packaging with FSC-certified kraft envelopes

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215745	HS Code	9018901000
Warranty (months)	12		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

### OTHER FEATURES

📅 08/06/2026 23:58:09 – For information purposes only, can be subject to change