

## REUSABLE NIBP FOR VITALIMATE MINI

Ref. MONI011A



CE



 09/06/2026 01:41:56 – For information purposes only, can be subject to change

# REUSABLE NIBP FOR VITALIMATE MINI

**Ref. MONI011A**

Vitalimate MINI reusable NIBP cuffs

## Description

The Vitalimate MINI NIBP cuffs are comfortable and come in all sizes for all types of patients.

- Reusable adult NIBP cuff (27-35cm) : ref MONI011
- Reusable large adult NIBP cuff (34-43cm) : ref MONI011A
- Small reusable adult NIBP cuff (20.5-28 cm): ref MONI011G
- Reusable pediatric NIBP cuff (14-21.5 cm): ref MONI011C
- Reusable neonatal NIBP cuff (10-15 cm): ref MONI011D
- Reusable adult NIBP thigh cuff (42-54 cm): ref MONI011B

Disposable neonatal cuffs are also available for infants.

## Features

CE Marking	CE	Medical device	Class I
Bar code	6948309002490	HS Code	9018193090
Warranty (months)	6	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.08
Bare Unit / Weight (kg)	0.075		

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

### OTHER FEATURES

📅 09/06/2026 01:41:56 – For information purposes only, can be subject to change