

SPENGLER ARMBAND, VAQUEZ-LAUBRY NANO POCKET

Ref. 330210



📅 10/06/2026 01:00:40 – For information purposes only, can be subject to change

SPENGLER ARMBAND, VAQUEZ-LAUBRY NANO POCKET



Ref. 330210

Description

Made of black Oxford-grade nylon, with inner bladder and dual tubing (10 cm & 60 cm), with split leather manometer holder, and Velcro fastening.

Child/Small adult (S), nylon, velcro, black : 330220

Adult (M), nylon, velcro, black : 330210

Adult (L), nylon, velcro, black : 330200

Features

CE Marking	NONE	Medical device	No
Bar code	3700446028596	HS Code	9018901000
Warranty (months)	6	Country	Indonesia

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.18

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 10/06/2026 01:00:40 – For information purposes only, can be subject to change