

WEIGHING PLATFORM WHEELCHAIR/SCALE MWA (300 KG) CLASS III VERIFIC

Ref. EKMWA01_CIII

KERN



CE
0297

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
3 ANS
3
3 YEARS
GUARANTEE

11/06/2026 01:09:20 – For information purposes only, can be subject to change

WEIGHING PLATFORM WHEELCHAIR/SCALE MWA (300 KG) CLASS III VERIFIC

KERN[®]

Ref. EKMWA01_CIII

Weighing platform for rolling wheelchair EKMWA01 approved in medical class III

Description

Weighing Platform wheelchair

- LCD display, digit height 25 mm
- Display Size (WxDxH) 200 x 55 x 128 mm
- Overall dimensions (LxWxH) 1180 x 1000 x 1130 mm
 - Weighing surface dimensions: L x D 880 x 840 mm
 - Weighing platform dimensions: L x D x H 995 x 945 x 40 mm
 - Column height: 940 mm
- Cable length of display approx. 180 cm
- Optional battery operation, 6 x 1.5V AA not included, operating time approx. 20 pm
- Power adapter (external) included • Net weight approx. 57.6 kg

 11/06/2026 01:09:20 – For information purposes only, can be subject to change

WEIGHING PLATFORM WHEELCHAIR/SCALE MWA (300 KG) CLASS III VERIFIC



Ref. [EKMWA01_CIII](#)

Features

CE Marking	CE0297	Medical device	Class Im (measuring function)
Bar code	3661040173627	HS Code	84231090
Warranty (months)	36	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	100
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	113
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	118
Packed unit / Weight	80	Bare Unit / Weight (kg)	68
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

11/06/2026 01:09:20 – For information purposes only, can be subject to change