

**WEIGHING PLATFORM WITH 2 BARRIERS (1 FOLDING BARRIER AND DISPLAY HOLDER AND 1 FIXED) MWA (300KG) CLASS III VERIFICATION**

**KERN®**

Ref. EKMWA03\_CIII



**CE**  
0297

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
**3**  
3 ANS  
3 YEARS  
GUARANTEE

11/06/2026 01:11:39 – For information purposes only, can be subject to change

# WEIGHING PLATFORM WITH 2 BARRIERS (1 FOLDING BARRIER AND DISPLAY HOLDER AND 1 FIXED) MWA (300KG) CLASS III VERIFICATION



[Ref. EKMWA03\\_CIII](#)

**Wheelchair weighing platform with 2 barriers (1 folding) and display door approved**

## Description

- Weighing Platform wheelchair
  - LCD display, digit height 25 mm
  - Display Size (WxDxH) 200 x 55 x 128 mm
  - Overall dimensions (LxWxH) 1180 x 1000 x 1130 mm
  - Weighing surface dimensions: L x D 880 x 840 mm
  - Weighing platform dimensions: L x D x H 995 x 945 x 40 mm
  - Column height: 940 mm
- Cable length of display approx. 180 cm
- Optional battery operation, 6 x 1.5V AA not included,  
operating time approx. 20 pm
- Power adapter (external) included • Net weight  
approx. 57.6 kg

 11/06/2026 01:11:39 – For information purposes only, can be subject to change

# WEIGHING PLATFORM WITH 2 BARRIERS (1 FOLDING BARRIER AND DISPLAY HOLDER AND 1 FIXED) MWA (300KG) CLASS III VERIFICATION



Ref. EKMWA03\_CIII

## Features

CE Marking	CE0297	Medical device	Class Im (measuring function)
Bar code	3661040173665	HS Code	84231090
Warranty (months)	36	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	100
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	113
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	118
Packed unit / Weight	94	Bare Unit / Weight (kg)	80
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 11/06/2026 01:11:39 – For information purposes only, can be subject to change