

Ref. 336401



CE



Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

12/06/2026 00:00:04 – For information purposes only, can be subject to change

=====

Ref. 336401



Features

| | | | |
|------------------------------------|---------------|---------------------------|------------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3700446026615 | HS Code | 9018901000 |
| Warranty (months) | 6 | Country | Indonesia |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.29 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

 12/06/2026 00:00:04 – For information purposes only, can be subject to change