

PULSE STETHOSCOPE / RED

Ref. 507522



CE



📅 12/06/2026 00:32:22 – For information purposes only, can be subject to change

PULSE STETHOSCOPE / RED

Ref. 507522



Description

For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaura with built-in spring.

Features

| | | | |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3700446009281 | HS Code | 90181990 |
| Warranty (months) | 36 | | |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.22 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Colors | Navy blue | Packed unit / Quantity | 1 |

📅 12/06/2026 00:32:22 – For information purposes only, can be subject to change