

Ref. BCTOCO34

COMEN



CE
1639

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE



14/06/2026 10:46:46 – For information purposes only, can be subject to change

=====

Ref. BCTOC034

COMeN

=====

Description

=====

Features

CE Marking	CE1639	Medical device	Class IIa
Bar code	6945429082470	HS Code	9018908400
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

 14/06/2026 10:46:46 – For information purposes only, can be subject to change