

LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

Ref. [MEEMER03_O](#)

EMERGENCY'S™



GARANTIE
2 ANS
2 YEARS
GUARANTEE



📅 14/06/2026 12:22:46 – For information purposes only, can be subject to change

LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

EMERGENCY'S™

Ref. [MEEMER03_O](#)

Intended to transport the most fragile + relief materials

Description

- Four compartments for interior supports and high-capacity compartment.
- Transparent pouch in the inner flap, designed to see and organize the material easily, thanks to the rubber-elastic that guarantees a correct fixation- Five outer pockets, two side pockets, one on the back and one on the top and a front pocket that can be divided internally thanks to rubber pockets and bands
- Base of the waterproof bag
- High resistance zips
- Reflective strips- Reinforced handles and shoulder straps
- External card holder
- Dimensions: 55 x 31.5 x 33 cm, capacity 38 L
- Internal dimensions: 40 x 20 x 24 cm
- Material 600D Polyester
- Weight 1.76 kg
- Sold empty

- Available in : Red (MEEMER03_R), Blue (MEEMER03_B) and Orange (MEEMER03_O)

 14/06/2026 12:22:46 – For information purposes only, can be subject to change

LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

EMERGENCY'S™

Ref. [MEEMER03_O](#)

Features

CE Marking	NONE	Medical device	No
Bar code	3661040202891	HS Code	420299000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	28	Bare Unit / Width (cm)	31.5
Packed unit / Height (cm)	20	Bare Unit / Height (cm)	33
Packed unit / Length (cm)	50	Bare Unit / Length (cm)	55
Packed unit / Weight	1.8	Bare Unit / Weight (kg)	1.76
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Orange	Packed unit / Quantity	1

📅 14/06/2026 12:22:46 – For information purposes only, can be subject to change