

**VAQUEZ-LAUBRY CLASSIC/VELCRO
CUFF/LARGE ADULT (L) / NAVY BLUE**

Ref. 518011



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
5 ANS
5
5 YEARS
GUARANTEE

📅 16/06/2026 01:15:23 – For information purposes only, can be subject to change

VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518011

=====

Description

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in nylon pouch.

Adult (M), cotton, velcro, navy: 518012

Adult (M), cotton, velcro, red: 518013

Adult (M), navy, with straps: 518112

Adult (M), nylon, velcro, grey: 518011

Adult (M), red with straps: 518113

Adult (L), nylon, velcro, grey: 518002

Child/small adult (S), cotton, velcro, motifs : 518022

VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518011

Features

| | | | |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking | CE0459 | Medical device | Class Im (measuring function) |
| Bar code | 3700446007287 | HS Code | 90189010 |
| Warranty (months) | 60 | Country | France |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 14 | Packed unit / Height (cm) | 8 |
| Packed unit / Length (cm) | 20 | Packed unit / Weight | 0.75 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 70 |
| Minimum humidity | 20 | Maximum humidity | 85 |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

16/06/2026 01:15:23 – For information purposes only, can be subject to change