

**F.O. MACINTOSH BLADE WITH INTEGRATED
LIGHT GUIDE, LUXAMED®,**

Ref. E1.413.012

LUXAMED®



17/06/2026 00:50:06 – For information purposes only, can be subject to change

F.O. MACINTOSH BLADE WITH INTEGRATED LIGHT GUIDE, LUXAMED[®],

LUXAMED[®]Ref. E1.413.012

LUXAMED MACINTOSH laryngoscope blades with integrated fibre optics.

Description

Laryngoscope blade with integrated optical glass

- High-quality stainless steel construction for a long service life
- No dust collecting cracks – easy cleaning and sterilisation
- Maintenance-free
- 8.000 fibres to achieve the maximum service life and illumination
- Optimal light emission due to the largest possible cross section (3 x 7 mm)
- Brightness approx. 14.000 Lux at 3.7 V LED and 2.5 V LED
- Compatible with all handles according to ISO 7376 (green standard)
- Autoclavable up to 134°C – approx. 4.000 times
- Data-Matrix with most important information

Differents size available :

0 : E1.410.012

1 : E1.411.012

2 : E1.412.012

3 : E1.413.012

4 : E1.414.012

5 : E1.415.012

**F.O. MACINTOSH BLADE WITH INTEGRATED
LIGHT GUIDE, LUXAMED®,**

LUXAMED®

Ref. E1.413.012

Features

| | | | |
|------------------------------------|----------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 04260381960522 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Germany |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | 5 | Max temperature | 40 |
| Minimum humidity | 10 | Maximum humidity | 90 |
| OTHER FEATURES | | | |
| Size | Size 3 | | |

17/06/2026 00:50:06 – For information purposes only, can be subject to change