


BABINSKY HAMMER, BLACK COLOR WITH PLASTIC HANDLE S.S

Ref. M05-320-280



17/06/2026 00:50:26 – For information purposes only, can be subject to change

BABINSKY HAMMER, BLACK COLOR WITH PLASTIC HANDLE S.S



Ref. M05-320-280



Description

- Handle with tip on the screw
- Unscrewable head 2 positions
- Length 25 cm
- Available in blue (M05-320-260), black (M05-320-250) and red (M05-320-280).

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040215471	HS Code	90189090
Warranty (months)	60	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	25
Packed unit / Weight	0.089	Bare Unit / Weight (kg)	0.089

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Red	Packed unit / Quantity	1
--------	-----	------------------------	---

📅 17/06/2026 00:50:26 – For information purposes only, can be subject to change