

OPERATING SCISSORS B/S STRAIGHT 11.5 CM

Ref. IC05411



OPERATING SCISSORS B/S STRAIGHT 11.5 CM



Ref. IC05411

Ideal for performing small treatments (dressings, compresses, etc.) without hurting the patient

Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance. The steels come into contact with human tissues, so they must comply with international standards and regulations. Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual

- Reference 100% IN-STORE: this service instrumentation reference is guaranteed to be available immediately.

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104741386	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Bare Unit / Width (cm)	4.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.4
Packed unit / Length (cm)	6	Bare Unit / Length (cm)	11.6
Packed unit / Weight	0.03	Bare Unit / Weight (kg)	0.036
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		