

KOCHER FORCEPS 1X2 TEETH STRAIGHT 16 CM

Ref. IP02016



KOCHER FORCEPS 1X2 TEETH STRAIGHT 16 CM



Ref. IP02016

Allows you to clamber or fix elements (tubing, etc.). Used mainly on sturdy fabrics.

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104742901	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Bare Unit / Width (cm)	7.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.4
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	16.5
Packed unit / Weight	0.05	Bare Unit / Weight (kg)	0.044
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		