

ANGULATED DISSECTING FORCEPS 13 CM (BRUCELLE)

Ref. IP04513



ANGULATED DISSECTING FORCEPS 13 CM (BRUCELLE)

Ref. IP04513



Allows to carry out care and facilitate the gripping of small objects

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104713529	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	18	Bare Unit / Width (cm)	0.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.6
Packed unit / Length (cm)	6	Bare Unit / Length (cm)	13
Packed unit / Weight	0.03	Bare Unit / Weight (kg)	0.016

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---