

## STERILIZING STAINLESS STEEL BOX 25X12X6 CM

Ref. AD25126



# STERILIZING STAINLESS STEEL BOX 25X12X6 CM

Ref. [AD25126](#)



**Convenient for storing surgical instruments**

## Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

## Features

|                   |                |                |          |
|-------------------|----------------|----------------|----------|
| CE Marking        | CE             | Medical device | Class I  |
| Bar code          | 08961104740044 | HS Code        | 90189084 |
| Warranty (months) | 12             | Country        | Pakistan |

### DIMENSIONS AND WEIGHT

|                           |       |                         |       |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm)  | 14    | Bare Unit / Width (cm)  | 12    |
| Packed unit / Height (cm) | 8     | Bare Unit / Height (cm) | 6     |
| Packed unit / Length (cm) | 27    | Bare Unit / Length (cm) | 25    |
| Packed unit / Weight      | 0.744 | Bare Unit / Weight (kg) | 0.726 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|