

**S/S KOCHER FORCEPS JAR 25 (DIAM)X H.
140 MM**

Ref. AI11615



S/S KOCHER FORCEPS JAR 25 (DIAM)X H. 140 MM

Ref. AI11615



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104740877	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	9	Bare Unit / Width (cm)	2.5
Packed unit / Height (cm)	9	Bare Unit / Height (cm)	14
Packed unit / Length (cm)	17	Bare Unit / Length (cm)	14
Packed unit / Weight	0.126	Bare Unit / Weight (kg)	0.121

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---