

S/S SERVING FORCEPS JAR WITH COVER(DIAM)50X200 MM

Ref. AI11723



19/06/2026 02:46:42 – For information purposes only, can be subject to change

S/S SERVING FORCEPS JAR WITH COVER(DIAM)50X200 MM

Ref. AI11723



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104740891	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	11	Bare Unit / Width (cm)	5
Packed unit / Height (cm)	11	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	24	Bare Unit / Length (cm)	20
Packed unit / Weight	0.222	Bare Unit / Weight (kg)	0.222
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		