

S/S KIDNEY DISH 24 X 13.5 X 4 CM

Ref. AI12544



19/06/2026 01:34:04 – For information purposes only, can be subject to change

S/S KIDNEY DISH 24 X 13.5 X 4 CM

Ref. AI12544



Ideal for storing instruments or temporarily installing medical waste

Description

- Decontamination with detergent and disinfectant -

Autoclavable (134°-18 minutes)- In the "

ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual

- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately

- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104739741	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	16	Bare Unit / Width (cm)	13.5
Packed unit / Height (cm)	10	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	28	Bare Unit / Length (cm)	25
Packed unit / Weight	0.187	Bare Unit / Weight (kg)	0.187

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---