

S. S BED PAN TYPE 37 X 29CM

Ref. AI13729



Description

Capacity 5 Liters

Maximum weight supported: 90 kg

19/06/2026 06:35:30 – For information purposes only, can be subject to change

S. S BED PAN TYPE 37 X 29CM

Ref. AI13729



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104739659	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	38	Bare Unit / Width (cm)	29
Packed unit / Height (cm)	32	Packed unit / Length (cm)	54
Bare Unit / Length (cm)	37	Packed unit / Weight	0.69
Bare Unit / Weight (kg)	0.68		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--