

S/S KIDNEY DISH 16 X 8.5 X 3.5CM

Ref. AI16124



19/06/2026 02:23:09 – For information purposes only, can be subject to change

S/S KIDNEY DISH 16 X 8.5 X 3.5CM

Ref. AI16124



Ideal for storing instruments or temporarily installing medical waste

Description

- Decontamination with detergent and disinfectant -

Autoclavable (134°-18 minutes)- In the "

ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual

- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately

- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104739727	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	17	Bare Unit / Width (cm)	8.5
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	3.5
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	16
Packed unit / Weight	0.077	Bare Unit / Weight (kg)	0.077

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---