

**PERFORATED BOXES, 18 X 8 X 4
CM, PERFORATED TOP / BOTTOM**

Ref. AI18084P



PERFORATED BOXES, 18 X 8 X 4 CM, PERFORATED TOP / BOTTOM

Ref. AI18084P



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104749351	HS Code	90189084
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	21	Bare Unit / Width (cm)	8
Packed unit / Height (cm)	10	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	23	Bare Unit / Length (cm)	18
Packed unit / Weight	0.248	Bare Unit / Weight (kg)	0.248

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---