

**S/S BOX PERFORATED 20 X 10 X 5CM
PERFORATED TOP / BOTTOM**

Ref. AI20105P



S/S BOX PERFORATED 20 X 10 X 5CM PERFORATED TOP / BOTTOM

Ref. AI20105P



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104749399	HS Code	90189084
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Bare Unit / Width (cm)	10.5
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	20.7
Packed unit / Weight	0.407	Bare Unit / Weight (kg)	0.407
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		