

S/S BOX 20 X 10 X 5 CM

Ref. AI20105



19/06/2026 07:55:21 – For information purposes only, can be subject to change

S/S BOX 20 X 10 X 5 CM

Ref. [AI20105](#)



Convenient for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104749382	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	11	Bare Unit / Width (cm)	10.6
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	5.1
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	20.5
Packed unit / Weight	0.397	Bare Unit / Weight (kg)	0.397

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---