

**S/S BOX PERFORATED 20 X 10 X 6 CM
PERFORATED TOP / BOTTOM**

Ref. AI20106P



S/S BOX PERFORATED 20 X 10 X 6 CM PERFORATED TOP / BOTTOM

Ref. [AI20106P](#)



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104749405	HS Code	90189084
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	6	Bare Unit / Width (cm)	6.9
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	3.2
Packed unit / Length (cm)	6	Bare Unit / Length (cm)	6.9
Packed unit / Weight	0.391	Bare Unit / Weight (kg)	0.391
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		