

STAINLESS STEEL BOX, 25 X 10 X 5 CM - PERFORATED TOP/BOTTOM

Ref. AI25105P



19/06/2026 07:07:03 – For information purposes only, can be subject to change

STAINLESS STEEL BOX, 25 X 10 X 5 CM - PERFORATED TOP/BOTTOM

Ref. [AI25105P](#)



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104740167	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	26	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	12	Bare Unit / Length (cm)	25
Packed unit / Weight	6.7		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---