

**S/S BOX PERFORATED 35 X 18 X 8 CM
PERFORATED TOP / BOTTOM**

Ref. AI35188P



S/S BOX PERFORATED 35 X 18 X 8 CM PERFORATED TOP / BOTTOM

Ref. AI35188P



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104740211	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	19	Bare Unit / Width (cm)	18
Packed unit / Height (cm)	9	Bare Unit / Height (cm)	8
Packed unit / Length (cm)	38	Bare Unit / Length (cm)	36
Packed unit / Weight	1.142	Bare Unit / Weight (kg)	1.142

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---