

**S/S BOX PERFORATED 40 X 20 X 9 CM
PERFORATED TOP / BOTTOM**

Ref. AI40209P



19/06/2026 06:41:47 – For information purposes only, can be subject to change

**S/S BOX PERFORATED 40 X 20 X 9 CM
PERFORATED TOP / BOTTOM****Ref. AI40209P****Features**

CE Marking	CE	Medical device	Class I
Bar code	08961104740228	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	21	Packed unit / Height (cm)	11
Packed unit / Length (cm)	43	Packed unit / Weight	1.51
Bare Unit / Weight (kg)	1.4		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---