

STAINLESS STEEL INSTRUMENT TRAY 40X30X4 CM

Ref. AI40304



18/06/2026 23:07:25 – For information purposes only, can be subject to change

STAINLESS STEEL INSTRUMENT TRAY 40X30X4 CM

Ref. AI40304



Ideal for storing small surgical instruments

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104739833	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	35	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	46	Bare Unit / Length (cm)	40
Packed unit / Weight	1	Bare Unit / Weight (kg)	0.9
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		