

INFANT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT2(13 TO 20 CM)



Ref. DS99506



Description

Infant cuff (13 to 20 cm)

Washable nylon cuff without pocket with Velcro fastening

19/06/2026 04:20:45 – For information purposes only, can be subject to change

INFANT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT2(13 TO 20 CM)



Ref. [DS99506](#)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040075259	HS Code	90189010
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	11.2
Packed unit / Weight	0.052	Bare Unit / Weight (kg)	0.052

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---