

# METAL TONGUE DEPRESSOR S/S CHILD CURVED 135° 13 CM - 29X22 MM

Ref. IA02400



## Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations. Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016). In the "

ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

19/06/2026 07:05:54 – For information purposes only, can be subject to change

# METAL TONGUE DEPRESSOR S/S CHILD CURVED 135° 13 CM - 29X22 MM



Ref. IA02400

## Features

CE Marking	CE	Medical device	Class I
Bar code	08961104730793	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15
Packed unit / Weight	0.03	Bare Unit / Weight (kg)	0.03
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

19/06/2026 07:05:54 – For information purposes only, can be subject to change