

# REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 3 MM), 15CM

Ref. [IA08409](#)



## Features

| CE Marking        | CE             | Medical device | Class I  |
|-------------------|----------------|----------------|----------|
| Bar code          | 08961104739055 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Pakistan |

### DIMENSIONS AND WEIGHT

|                           |       |                         |       |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm)  | nc    | Bare Unit / Width (cm)  | 18    |
| Packed unit / Height (cm) | nc    | Bare Unit / Height (cm) | 0.2   |
| Packed unit / Length (cm) | nc    | Bare Unit / Length (cm) | 19    |
| Packed unit / Weight      | 0.015 | Bare Unit / Weight (kg) | 0.015 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|