


**BABINSKI REFLEX HAMMER, ADULT, 25 CM,
METAL**

Ref. ID0660_G



19/06/2026 05:19:45 – For information purposes only, can be subject to change

BABINSKI REFLEX HAMMER, ADULT, 25 CM, METAL



[Ref. ID0660_G](#)

Convenient for testing tendon reflexes

Description

- Handle with point on screw
- 2-position removable head
- Length: 25 cm

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104715622	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	28	Bare Unit / Width (cm)	4.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4.4
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	25
Packed unit / Weight	0.15	Bare Unit / Weight (kg)	0.08

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 19/06/2026 05:19:45 – For information purposes only, can be subject to change