

MAYO-HEGAR NEEDLE HOLDER 14 CM

Ref. IP29414



MAYO-HEGAR NEEDLE HOLDER 14 CM

Ref. **IP29414**



Allows you to hold a needle firmly during wound suturing

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104744370	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	23	Bare Unit / Width (cm)	9
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	14
Packed unit / Weight	0.04	Bare Unit / Weight (kg)	0.042

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---