

## FEICHENFELD SPLINT FORCEPS 10 CM

Ref. IP41410



19/06/2026 07:34:39 – For information purposes only, can be subject to change

# FEICHENFELD SPLINT FORCEPS 10 CM

Ref. IP41410



**Used during surgical procedures or care. Wide part to facilitate the handling.**

## Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

## Features

CE Marking	NONE	Medical device	No
Bar code	08961104714335	HS Code	90189084
Warranty (months)	24	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	18	Bare Unit / Width (cm)	1.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.1
Packed unit / Length (cm)	6	Bare Unit / Length (cm)	10
Packed unit / Weight	0.03	Bare Unit / Weight (kg)	0.018

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 19/06/2026 07:34:39 – For information purposes only, can be subject to change