

FOERSTER FORCEPS STRAIGHT 24 CM

Ref. IP43424



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

20/06/2026 00:57:22 – For information purposes only, can be subject to change

FOERSTER FORCEPS STRAIGHT 24 CM

Ref. IP43424



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104720138	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	24
Packed unit / Weight	0.078	Bare Unit / Weight (kg)	0.077
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 20/06/2026 00:57:22 – For information purposes only, can be subject to change