

COMPLETE NORMAL DELIVERY SET

Ref. IBOI105



18/06/2026 22:25:57 – For information purposes only, can be subject to change

COMPLETE NORMAL DELIVERY SET

Ref. IBOI105



Description

This box contains:

- Stainless steel box, 40 x 20 x 9 cm
- Tarnier forceps, with retractor, 40 cm
- Mousse scissors, curved, 18 cm, Holtex
- Jean-Louis Faure forceps, curved, A/G, 22 cm
- Pinard obstetric stethoscope, aluminum
- Bladder catheter, metal, female, 15 cm
- Bar umbilical clamp, 8 cm
- Membrane piercer, 19 cm
- Ribemont insufflator, without bulb, 10 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104746923	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.2
Bare Unit / Weight (kg)	3.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

18/06/2026 22:25:57 – For information purposes only, can be subject to change