

TRACHEOSTOMY +CANNULAS SURGERY SET

Ref. IBOI360



TRACHEOSTOMY +CANNULAS SURGERY SET

Ref. IBOI360



Description

This box contains:

- Laborde retractor, 3 prongs, for ENT
- Tracheal retractor, foam, 2 teeth, 16 cm (Joseph)
- Stainless steel box, 20 x 10 x 5 cm -IS
- Dissection forceps, S/G, 14 cm -IS, Holtex
- Tracheotome, 7 mm

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104748293	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.65
Bare Unit / Weight (kg)	0.65		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---