

RIB SHEAR 20 CM DOYEN TYPERef. IC32420**Features**

CE Marking	CE	Medical device	Class Ir
Bar code	08961104731493	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	21
Packed unit / Weight	0.253	Bare Unit / Weight (kg)	0.253

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---