

# HOOK TO REMOVE IUD 25CM (DOWNS P-347)

Ref. [IC38201](#)



## Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104738034	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.026
Bare Unit / Weight (kg)	0.024		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

19/06/2026 00:22:45 – For information purposes only, can be subject to change