

GINESTET HOOKRef. **IC39110****Features**

CE Marking	CE	Medical device	Class Ir
Bar code	08961104722163	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	2
Packed unit / Length (cm)	24	Packed unit / Weight	0.13
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		