


CHIROPODIST CURETTE FENESTRED 3 MM

Ref. IC64003



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

📅 19/06/2026 01:46:14 – For information purposes only, can be subject to change

CHIROPODIST CURETTE FENESTRED 3 MM

Ref. IC64003



For small pedicure procedures

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

Features

| | | | |
|-------------------|----------------|----------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 08961104725904 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm) | nc | Bare Unit / Width (cm) | 0.7 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 0.3 |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 13.5 |
| Packed unit / Weight | 0.018 | Bare Unit / Weight (kg) | 0.017 |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

📅 19/06/2026 01:46:14 – For information purposes only, can be subject to change