

CORWIN TONSIL FORCEPS

Ref. IC99300



CORWIN TONSIL FORCEPS

Ref. IC99300



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104730915	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16.6
Packed unit / Weight	0.083	Bare Unit / Weight (kg)	0.083

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---