

DOYEN MOUTH-GAG CHILD, 12 CM

Ref. ID03600



DOYEN MOUTH-GAG CHILD, 12 CM

Ref. ID03600



Description

Doyen children's mouth opener in stainless steel with a total length of 12 cm (Hooks: 45 mm - Ends: 21 mm).

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104730663	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	12
Packed unit / Weight	0.081	Bare Unit / Weight (kg)	0.081
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		