

DOYEN MOUTH-GAG ADULT, 13.5 CM

Ref. ID03602



DOYEN MOUTH-GAG ADULT, 13.5 CM

Ref. ID03602



Description

Dean stainless steel mouth opener for adults. Total length 13.5 cm (Hooks: 60 mm - Ends: 24 mm).

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104730670	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14.3
Packed unit / Weight	0.121	Bare Unit / Weight (kg)	0.121

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---