

THOMAS WAX MODELING INSTRUMENTS N°3

Ref. IDINTH03



CE



18/06/2026 23:32:54 – For information purposes only, can be subject to change

THOMAS WAX MODELING INSTRUMENTS N°3



Ref. IDINTH03

Features

| | | | |
|------------------------------------|----------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class Ir |
| Bar code | 08961104739482 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 1 | Bare Unit / Width (cm) | 1 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 0.8 |
| Packed unit / Length (cm) | 18 | Bare Unit / Length (cm) | 18 |
| Packed unit / Weight | 0.028 | Bare Unit / Weight (kg) | 0.028 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

18/06/2026 23:32:54 – For information purposes only, can be subject to change